

Development of Novel Carrier Systems with Controlled Surface Features for Vascular Delivery

Joerg Lahann

Department of Chemical Engineering, University of Michigan, Ann Arbor,
2300 Hayward St., Ann Arbor, MI 48109, USA

Institute of Functional Interfaces, Karlsruhe Institute of Technology,
Hermann-von-Helmholtz-Platz 1, Eggenstein-Leopoldshafen, 76344, Germany
lahann@umich.edu

Extended Abstract

Nano- and microparticles have been successfully utilized in a variety of research areas, such as diagnostics, drug delivery or regenerative medicine. Compositional flexibility, tunability, and versatility of polymer particles has been accomplished by progress in particle fabrication methods. In particular, the complementary control of internal (bulk) and external (surface) features has been increasingly recognized as important design parameters for multifunctional particles.

In the past, particles have been fabricated with several distinct methods. In particular, electrospraying of polymer solutions is a straightforward method to fabricate nano-/microparticles as well as fibers by applying high electrical voltage to polymeric solutions. Building on these processes, electrohydrodynamic (EHD) co-jetting involves two or more capillary needles in a side-by-side configuration that allow different polymer solutions to be processed in parallel. Under the laminar flow regimen encountered in these systems, biphasic droplets are formed at the outlet point of the adjacent needles. Application of an electric field to the nozzle leads to a stable Taylor cone at the tip of the biphasic droplet.

Initial work has been focused on water-soluble polymers because environmentally friendly water-based jetting systems typically exhibit lower toxicity, which are important for biomedical applications. More recently, EHD co-jetting of organic-soluble polymers have been used to fabricate anisotropic microstructures. In particular, the EHD co-jetting of organic solution of poly(lactic-co-glycolic acid) (PLGA) has been successfully used to prepare biodegradable bicompartamental particles and fibers. Because of their biodegradability, these particles may constitute promising building blocks for bio-hybrid materials or find potential applications in drug delivery and tissue engineering. The systematic investigation of EHD co-jetting of PLGA reveals that the jetting solution and process parameters can control compartmentalization, shape, and size.

Multicompartamental particles feature dissimilar materials and combine orthogonal sets of properties in the same particle. Selective surface modification of parts of the surface can create surface patterns. The precise nano- and microscale control of the architecture of biodegradable materials is required for many biomedical applications, including controlled drug delivery for vascular applications.